

CELL CULTURE

From Sample Collection to LC-MS Analysis of Moderately Polar Intracellular Metabolites

STEPS 1-2




SAMPLE COLLECTION

1. Grow cells with 5-6 replicates/condition to yield ~1 mg protein/sample
2. Wash cells twice with room temperature PBS

Note: Remove PBS thoroughly as salts are not MS-friendly


STEPS 3-8



EXTRACTION

3. Place samples on dry ice and add cold MeOH:H₂O (4:1) to wells; scrape cells and transfer to a tube with beads
4. Vortex for 30 s
5. Homogenize
6. Incubate for 1 h at 20°C; centrifuge for 15 min at 4°C
7. Transfer supernatants to a new tube and evaporate solvent
8. Quantify protein concentration of cell pellets

STEPS 9-11



ANALYSIS

9. Reconstitute extracts in MeOH:H₂O (1:1)
10. Sonicate for 1 min in ice-cold bath; centrifuge for 15 min at 4°C
11. Transfer supernatant to UPLC vials and inject for LC-MS analysis

SAMPLE COLLECTION

- ① Grow cells in individual 100 or 60 mm plates, 5-6 biological replicates per sample group, to have an equivalent of ~1 mg of proteins or 1×10^6 cells per sample (*an estimation that depends on the type/size of cells*) → protocol can be adjusted for metabolite profiling from lower amounts of biological material (down to 100 K cells)
- ② Washing off the spent media: Wash the cells twice with PBS at room temperature (**DO NOT** expose your cells to cold stress using ice-cold solvent!) ~4 ml for 60 mm plates, ~10 ml for 100 mm plates (5-10s), to wash off the spent media x 2. Amino acids (and other metabolites) are usually abundant in the media and can contribute to the measurement of intracellular concentrations. Remove the PBS thoroughly (salts are NOT MS-friendly).

METABOLITE EXTRACTION

- ③ Place the cell cultures on dry ice and add 1 ml of cold MeOH:H₂O (4:1, v/v) (UPLC-MS grade) to quench the metabolism and burst/release the cell contents. Scrape the cells off the dish with cell scraper and transfer the solution with cells to a lysis tube and add the ceramic beads.
- ④ Vortex briefly (~30 seconds) and sonicate on ice at least 30 seconds to extract metabolites out of cells. This step is also required if performing freeze-thaw.
- ⑤ Homogenize the solution in tissue homogenizer (Precellys®) 3 x 20s at 6,000 rpm (with 5s between each homogenization). A cold trap (Cryolys) of the homogenizer should be filled with dry ice to keep the sample cold. **Remark:** If you don't have a tissue homogenizer you can perform 2 x freeze-thaw cycles to extract the cellular content. This is particularly important for bacterial cells.
- ⑥ Incubate samples for 1 hour at -20°C, then centrifuge at 16,000 g at 4°C for 15 minutes
- ⑦ Transfer the supernatants to another 1.5 ml Eppendorf and evaporate the supernatant to dryness in a vacuum concentrator (Labconco™ CentriVap Benchtop) at 10°C (*If the time is limited, the dry extracts can be stored at -80°C at this stage, prior to further processing*). Save cell pellets for step 8.

- ⑧ Quantify the total protein content on cell pellets (using BCA or another protein assay)
- ⑨ Reconstitute the dry extracts in MeOH/water (1:1, v/v) normalized to protein content
- ⑩ Sonicate for 1 minute in an ice-cold bath and centrifuge for 15 minutes at 16,000 g and 4°C to remove debris
- ⑪ Transfer supernatants to UPLC vials with inserts and inject them for LC/MS analysis. (The storage at -80°C should be avoided at this stage because it might cause the metabolite precipitation and therefore compromise the coverage)

LIQUID CHROMATOGRAPHY CONDITIONS

- LC System: UPLC
- Column & Guard: WATER ACQUITY UPLC BEH C18 Column, 130Å, 1.7 µm, 2.1 mm X 100 mm BEH C18, 2.1x 100 mm, 1.7 µm & WATER ACQUITY UPLC BEH C18 VanGuard Pre-column, 130Å, 1.7 µm, 2.1 mm X 5 mm
- Flow rate: 400 ul/min
- Mobile phases:

| Time (minutes) | A (H ₂ O + 0.1% FA) | B (ACN + 0.1% FA) |
|----------------|--------------------------------|-------------------|
| 0 | 99 | 1 |
| 1 | 99 | 1 |
| 10 | 1 | 99 |
| 13 | 1 | 99 |
| 14 | 99 | 1 |
| 17 | 99 | 1 |

MASS SPECTROMETRY CONDITIONS

- Ionization modes: ESI (+) and ESI (-)
- Acquisition mode: MS and Auto MS/MS
- m/z range: 50-1,000 Da
- Spectra rate: 500 milliseconds per MS spectra, 130 milliseconds per MS/MS spectra (10-50 eV)

MATERIALS

- All solvents are HPLC grade
- UPLC Column Waters ACQUITY PREMIER BEH C18, 130Å, 1.7 µm VanGuard FIT 2.1 x 100 mm MVK
- Glass vials Thermo Fisher, Part No. 4000-S1W
- Vial caps Agilent, Part No. 5182-0717
- Glass inserts Agilent, Part No. 5181-8872
- 1.5 ml tubes (natural) Eppendorf, Part No. 022363204
- Standards mixture Cayman, Vascular Eicosanoid Urinary Metabolite LC-MS Mixture, Part No. 19668

REFERENCE

<https://www.nature.com/articles/s41592-021-01197-1>
<https://www.nature.com/articles/s41592-020-0942-5>

RESOURCES

METLIN Gen2 (<https://massconsortium.com/>)
 XCMS Online (<https://xcmsonline.scripps.edu/>)
 Cayman Chemical Lipidomic Standards (<https://www.caymanchem.com/>)